

BOSTON SCOLIOSIS BRACE™

Built-in iliac crest pads allow tighter pressure of the intervertebral spine. A vertical anterior center line is marked on all unfinished braces. Module design incorporates symmetrical modification and uniform dimensions according to size. Standard size modules will fit 80% of all adolescent scoliosis cases up to level T-6.

Size	Hips cm	Waist cm	Xyphoid cm	ASIS cm	Model No.	
Juvenile: Std. 5/32" Copolymer					0° Lordosis	15° Lordosis
J1	54	40	49	14	00J1	15J1
J2	58	42	51	15	00J2	15J2
J3	62	45	54	16	00J3	15J3
J4	64	47	56	17	00J4	15J4
Youth: Std. 11/64" Copolymer						
Y1A	66	51	55	17	00Y1A	15Y1A
Y1	67	46	56	18	00Y1	15Y1
Y2	71	47	60	19	00Y2	15Y2
Y2A	72	55	62	20	00Y2A	15Y2A
Y3	75	48	63	20	00Y3	15Y3
Y3A	78	57	69	21	00Y3A	15Y3A
Y4	79	51	64	23	00Y4	15Y4
Y4F	80	57	66	22	00Y4F	15Y4F
Y5	82	52	68	24	00Y5	15Y5
Y5F	82	58	68	23	00Y5F	15Y5F
Y6	84	56	73	24	00Y6	15Y6
Y6A	84	63	78	23	00Y6A	15Y6A
Y7	86	54	69	25	00Y7	15Y7
Y7F	86	60	69	24	00Y7F	15Y7F
Y7A	86	69	73	24	00Y7A	15Y7A
Y8	88	58	75	25	00Y8	15Y8
Y8F	88	64	75	24	00Y8F	15Y8F
Y8A	90	72	83	24	00Y8A	15Y8A
Y9	92	60	72	26	00Y9	15Y9
Y9A	92	63	69	25	00Y9A	15Y9A
Y9F	92	68	78	25	00Y9F	15Y9F
Adolescent: Std. 3/16" Copolymer						
A1	94	62	76	27	00A1	15A1
A2	94	70	84	26	00A2	15A2
A2A	96	66	75	25	00A2A	15A2A
A3	97	68	80	27	00A3	15A3
A4	97	76	88	27	00A4	15A4

* Modules with a hip circumference less than 54 cm will be made of 1/8" copolymer.